



COMPUTER VISION

Presented by Anna D'Angela | December 15th, 2020

Phase 4 Project - Neural Networks


A decorative vertical bar on the left side of the slide, transitioning from purple at the top to blue at the bottom. In the top right corner, there are several concentric circular lines in a light blue color.

BUSINESS CASE

Even for a trained radiologist, it is challenging to examine chest X-rays for markers of pneumonia. There is a need to improve the diagnostic accuracy.

A medical imaging company is looking to increase efficiency for its technicians by building a classifier to aid them in detecting pneumonia from chest X-ray images.

Competition Data Hosted by Kaggle
Project for Flatiron School, Online Data Science Program




X-RAY CLASSIFICATION OF PEDIATRIC PNEUMONIA

“Pediatric pneumonia is responsible for the deaths of more than **800,000** young children worldwide each year.”


—United Nations Children's Fund
(UNICEF)

EXAMPLE X-RAYS


Normal Lungs, 1




Normal Lungs, 2



Normal Lungs, 3



Normal Lungs, 4




Normal Lungs, 5



Pnuemonia Lungs, 1




Pnuemonia Lungs, 2




Pnuemonia Lungs, 3



Pnuemonia Lungs, 4



Pnuemonia Lungs, 5



CLINICAL EFFICIENCY



PATIENT VISIT

Physical exam and chest x-ray if symptoms present



RAPID ANALYSIS

Radiologist is assisted in diagnosis by the model



TREATMENT

A swift start to treatment leads to faster recovery

The background is a solid dark purple color. On the left side, there are several concentric white circles of varying radii. A single horizontal white line spans the width of the slide, positioned above the 'OBJECTIVE' section.

OBJECTIVE

Increase clinical efficiency in diagnosis of pediatric pneumonia

METHOD

Create an X-ray image classification tool to assist in labeling cases of pediatric pneumonia

SUCCESS CRITERIA

Optimize for accuracy to limit false predictions

89.77%

MODEL DIAGNOSTIC ACCURACY

90%

Accuracy

94%

Recall

82%


Precision

0.90

F1 Score


MODEL PERFORMANCE

PREDICTIONS




‘NORMAL’

PREDICTIONS




‘NORMAL’

PREDICTIONS



‘PNEUMONIA’

PREDICTIONS



'PNEUMONIA'

- Continue collecting labeled images to progressively train the model.
- Store image data at 128 x 128 to conserve storage memory (this is up to a 10% reduction in original image size).
- Use the model to improve efficiency of X-ray review, but do not replace human classification.

RECOMMENDATIONS

FUTURE WORK

MODEL TUNING

- Increase the quantity of images in the training set.
- Utilize transfer learning to improve the base model.
- User experience: Build an application to receive an X-ray as input and output a prediction.

THANK YOU!

Please find my full analysis on
GitHub: @anna-dang

Anna D'Angela | Detroit, MI

CREDITS: This presentation template was
created by Slidesgo, including icons by Flaticon,
and infographics & images by Freepik.

APPENDIX (i)

CLASSIFICATION REPORT

Classification Report for Test Data:

	precision	recall	f1-score	support
0.0	0.90	0.82	0.86	176
1.0	0.90	0.94	0.92	293
accuracy			0.90	469
macro avg	0.90	0.88	0.89	469
weighted avg	0.90	0.90	0.90	469


15/15 [=====] - 0s 7ms/step - loss:
Loss of the model is - 0.2689375579357147

15/15 [=====] - 0s 7ms/step - loss:
Accuracy of the model is - 89.76545929908752 %

Correct: 421, 89.77%
Incorrect: 48, 10.23%

ROC CURVE / AUC

A perfect binary classifier has an 'Area Under Curve' (AUC) of 1, with ROC curve hugging the top left corner



APPENDIX (ii)

IMAGE PREPROCESSING

- Images were resized to 124 x 124 pixels, with 3 RGB color channels
- Pixel values were normalized to a 0-1 scale
- To prepare the model to discern noise, four data augmentations were used: rotation, vertical/horizontal shifting, and zoom
- The imbalanced data set (75% pneumonia vs. 25% normal X-rays) was corrected by applying class weights


MODEL STRUCTURE

Model: "sequential_15"

Layer (type)	Output Shape	Param #
conv2d_64 (Conv2D)	(None, 128, 128, 128)	3584
dropout_8 (Dropout)	(None, 128, 128, 128)	0
max_pooling2d_55 (MaxPooling)	(None, 64, 64, 128)	0
conv2d_65 (Conv2D)	(None, 62, 62, 64)	73792
dropout_9 (Dropout)	(None, 62, 62, 64)	0
max_pooling2d_56 (MaxPooling)	(None, 31, 31, 64)	0
conv2d_66 (Conv2D)	(None, 29, 29, 32)	18464
dropout_10 (Dropout)	(None, 29, 29, 32)	0
max_pooling2d_57 (MaxPooling)	(None, 14, 14, 32)	0
flatten_15 (Flatten)	(None, 6272)	0
dense_30 (Dense)	(None, 32)	200736
dense_31 (Dense)	(None, 1)	33
Total params: 296,609		
Trainable params: 296,609		
Non-trainable params: 0		



APPENDIX (iii)

3D CNN MODEL ARCHITECTURE



APPENDIX (iv)

IMAGE EXPLAINER (LIME)



APPENDIX (iv)

IMAGE EXPLAINER (LIME)

