



# COMPUTER VISION

Presented by Anna D'Angela | December 15th, 2020

Phase 4 Project - Neural Networks


A decorative vertical bar on the left side of the slide, transitioning from purple at the top to blue at the bottom. In the top right corner, there are several concentric, dashed blue circles of varying sizes.

## BUSINESS CASE

Even for a trained radiologist, it is challenging to examine chest X-rays for markers of pneumonia. There is a need to improve the diagnostic accuracy.

A medical imaging company is looking to increase efficiency for its technicians by building a classifier to aid them in detecting pneumonia from chest X-ray images.

Competition Data Hosted by Kaggle  
Project for Flatiron School, Online Data Science Program



# **X-RAY CLASSIFICATION OF PEDIATRIC PNEUMONIA**

“Pediatric pneumonia is  
responsible for the deaths  
of more than **800,000**  
young children worldwide  
each year.”

—United Nations Children's Fund  
(UNICEF)

# EXAMPLE X-RAYS

Normal Lungs




Normal Lungs



Normal Lungs



Normal Lungs



Normal Lungs



Pnuemonia Lungs



Pnuemonia Lungs




Pnuemonia Lungs



Pnuemonia Lungs



Pnuemonia Lungs



# CLINICAL EFFICIENCY

---



## PATIENT VISIT

Physical exam and chest x-ray if symptoms present



## RAPID ANALYSIS

Radiologist is assisted in diagnosis by the model



## TREATMENT

A swift start to treatment leads to faster recovery

---

## OBJECTIVE

Create an X-ray image classification tool to assist in diagnosing cases of pediatric pneumonia


## METHOD

Use a convolutional neural network (CNN) to find patterns in X-rays to predict class label

## SUCCESS CRITERIA

Optimize for accuracy to limit false predictions

# CNN MODEL ARCHITECTURE





# 89.77%

---

**MODEL DIAGNOSTIC ACCURACY**

90%

Accuracy

94%

Recall

82%


Precision

0.90

F1 Score


**MODEL PERFORMANCE**

# PREDICTIONS




**‘NORMAL’ /  
NEGATIVE CASES**

# PREDICTIONS




**'NORMAL' /  
NEGATIVE CASES**

# PREDICTIONS



**‘PNEUMONIA’ /  
POSITIVE CASES**

# PREDICTIONS



**‘PNEUMONIA’ /  
POSITIVE CASES**

# USER RECOMMENDATIONS

- **Model use:** Use the model as a pre-screening tool to improve efficiency of X-ray review, but do not replace human classification.
- Continue collecting labeled images to progressively train the model.
- Store image data at 128 x 128 to conserve storage memory (this is up to a 10% reduction in original image size).

# FUTURE WORK

## MODEL TUNING / IMPROVEMENTS

---

- Increase the quantity of images in the training set.
- Utilize transfer learning to improve the base model.
- **User experience:** Build an application to receive an X-ray as input and output a prediction.



# THANK YOU!

Please find my full analysis on  
GitHub: @anna-dang

Anna D'Angela | Detroit, MI

CREDITS: This presentation template was  
created by Slidesgo, including icons by Flaticon,  
and infographics & images by Freepik.

# APPENDIX (i)

## CLASSIFICATION REPORT

-----  
Classification Report for Test Data:

	precision	recall	f1-score	support
0.0	0.90	0.82	0.86	176
1.0	0.90	0.94	0.92	293
accuracy			0.90	469
macro avg	0.90	0.88	0.89	469
weighted avg	0.90	0.90	0.90	469


-----  
15/15 [=====] - 0s 7ms/step - loss:  
Loss of the model is - 0.2689375579357147

15/15 [=====] - 0s 7ms/step - loss:  
Accuracy of the model is - 89.76545929908752 %

-----  
Correct: 421, 89.77%  
Incorrect: 48, 10.23%  
-----

## ROC CURVE / AUC

A perfect binary classifier has an 'Area Under Curve' (AUC) of 1, with ROC curve hugging the top left corner



# APPENDIX (ii)

## IMAGE PREPROCESSING

- Images were resized to 124 x 124 pixels, with 3 RGB color channels
- Pixel values were normalized to a 0-1 scale
- To prepare the model to discern noise, four data augmentations were used: rotation, vertical/horizontal shifting, and zoom
- The imbalanced data set (75% pneumonia vs. 25% normal X-rays) was corrected by applying class weights


## MODEL STRUCTURE

Model: "sequential\_15"


Layer (type)	Output Shape	Param #
conv2d_64 (Conv2D)	(None, 128, 128, 128)	3584
dropout_8 (Dropout)	(None, 128, 128, 128)	0
max_pooling2d_55 (MaxPooling)	(None, 64, 64, 128)	0
conv2d_65 (Conv2D)	(None, 62, 62, 64)	73792
dropout_9 (Dropout)	(None, 62, 62, 64)	0
max_pooling2d_56 (MaxPooling)	(None, 31, 31, 64)	0
conv2d_66 (Conv2D)	(None, 29, 29, 32)	18464
dropout_10 (Dropout)	(None, 29, 29, 32)	0
max_pooling2d_57 (MaxPooling)	(None, 14, 14, 32)	0
flatten_15 (Flatten)	(None, 6272)	0
dense_30 (Dense)	(None, 32)	200736
dense_31 (Dense)	(None, 1)	33
Total params: 296,609		
Trainable params: 296,609		
Non-trainable params: 0		

## APPENDIX (iii)

### IMAGE EXPLAINER (LIME)



**'NORMAL'**



## APPENDIX (iii)

### IMAGE EXPLAINER (LIME)

